



Oasis Sanic Reverse Osmosis System

Antimicrobial Water Purification System with Bacteria Static Technology
For Drinking Water, Ice, Hot Water, and Cooking Applications

Atlas Filtri® Sanic, Under Counter Reverse Osmosis Water Purification System is designed to meet all your drinking and cooking water requirements.

The system is designed to provide up to 50 gallons of pure water a day. In addition to the secondary faucet at the sink, the system can be connected to a refrigerator for ice and water dispensing, an under counter ice machine, a hot water dispensing system for all four points of use.

The system will reduce dissolved contaminants in the water supply up to 99% through the RO membrane separation process. The ultra-fine membrane filters everything larger than 0.0001-micron, filtering everything except pure water.

Examples of contaminants removed are chlorine, fluoride, arsenic, aluminum, barium, cadmium, chlorides, chromium, lead, mercury, and viruses. The patented single-molded manifolds with three housings reduce potential leak points and make the system more efficient.



The clear sediment first stage filter gives an indication of just how dirty the water can get over months of water production. Antimicrobial water contact components are designed to eliminate the growth of bacteria that may be present in the water supply.

Water after an RO membrane can have a low pH level when certain minerals are removed. The alkaline re-mineralization cartridge raises the pH of the water to improve the taste.

Atlas Filtri® is driven to providing our customer with the best available products to meet specific filtration requirements. This is done by using the most advanced manufacturing and design methods. The international patents received come from a constant commitment to research and develop that result in new and innovative products.

Discover more at us.atlasfiltri.com.



TW30-1812-50HR elements have been tested and certified by **NSF** International against **NSF/ANSI** Standard 58 for material requirements



Housings certified by IAPMO RWT to NSF/ANSI 42, 61, 372 lead free, CSA B483.1. In compliance with drinking water regulations, USA and Canada

Replacement Parts for Reverse Osmosis System

Filter change-outs vary depending upon type



Oasis Sanic Pre-Filter Cartridge Kit

Stage 1: Sediment Antimicrobial 25-micron filter
Stage 2: Sediment Antimicrobial 5-micron filter
Stage 3: VOC/CYST 0.5-micron carbon block filter
Change Out Every 12 months

O-ring Kit

Three O-rings replaced on filter change-out
Order with Pre-Filter Cartridge Kit



LubriKit

O-ring lubricant required for proper seal
Order with Pre-Filter Cartridge Kit



FilterFIT

Pre-Filter Centering Device for proper alignment
Order with Pre-Filter Cartridge Kit



Pre-Filter Housing Wrench

Device Designed to Open Pre-Filter Housing for Filter Cartridge Replacement
Available Accessory




Post-Inline Filter

Activated Carbon Polishing Filters
Change Out Every 12 months



Remineralization Alkaline Cartridge

Adjusts the pH of the water to improve taste
Change Out Every 12 months



Reverse Osmosis Membrane

Filter Cartridge – 50 gallons per day
Change Out every 36 months

Part Number	Description	Barcode
RE7402010	RO Membrane 50 GPD – 1812S	8021647229557
RE7402004	Post Inline Activated Carbon Filter	8021647229496
RE7402005	Alkaline Re-Mineralization Cartridge	8021647229502
SE6075221	Sanic Pre-Filter Cartridge Kit (Quantity 3)	8021647284525
RE7514331	O-ring Kit for DP TRIO (Quantity 3)	
RE9020660	LUBRIKit O-ring Lubricant	8021647270924
LB7406012	FilterFIT Filter Cartridge Centering Device	8021647137982
RB7403017	Wrench RO System DP	8021647164094

Check out more about Atlas Filtri at us.atlasfiltri.com

Atlas Filtri® is a registered trademark of ATLAS FILTRI srl. Unauthorized use of the registered trademark is prohibited. Images and context are the property of ATLAS FILTRI srl, which reserves the right to change product design and specification without prior notification.

Atlas Filtri North America LLC

1068 North Farms Road, Building 3, Wallingford, CT 06492, USA

Office +1 (203) 284-0080 - Fax +1 (203) 294-9226

email: atlasfiltrinorthamerica@atlasfiltri.com

A Company of the Atlas Filtri Group. World headquarters and production facilities located in Italy.

